TERMOMETRO FRONTALE AD INFRAROSSI **THERMO-PRECISION**

Infrared forehead thermometer

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Thermo-Precision infrared forehead thermometer

Dear customer,

We thank you for having purchased a Thermo-Precision Forehead, the thermometer that takes a fast, precise body temperature reading from a distance.

Thermo-Precision Forehead is an innovative method for measuring body temperature. Thanks to the exclusive infrared technology it detects body temperature from a distance of 0.5 cm (0.2 inch) up to the distance of 10 cm (4 inch) from the forehead.

Innovative technology and modern design of Thermo-Precision Forehead is the solution, ideal for the entire family.

The thermometer measures the infrared radiation given off by the surface of the forehead and converts this measurement into the corresponding oral temperature.

Thermo-Precision Forehead is produced using the most advanced technologies and undergo rigorous quality controls, in compliance with EC regulations, thus guaranteeing utmost quality and precision.

BODY TEMPERATURE

Body temperature varies from person to person and undergoes different fluctuations during the course of the day, running approximately from 35.5°C to 37.8°C (95.5°F - 100.4°F).

For these reasons, to correctly determine the temperature, it is best to know one's normal, healthy temperature as measured on the forehead at various times of the day. This will provide a valid reference when evaluating any real rise in temperature. To determine if one has a fever, compare the temperature detected with the person's normal temperature. Remember, a rise over the reference body temperature of 1° C (1.8°F) or more is generally indication of fever.

HOW TO TAKE MEASUREMENT

Tips for a correct temperature measurement

In order to take the temperature in the best possible way, a few simple rules should be followed:

- The thermometer automatically adjusts its reading of forehead temperature on the basis of room temperature. For this reason, in order to obtain a correct reading, the thermometer should be in the room where measurement will take place at least 20-25 minutes prior to measurement. Also the person, whose temperature will be taken, should remain in the same room as the thermometer for at least 15 minutes before having his/her temperature taken.
- Before taking the temperature, wipe the forehead to remove perspiration and wait a few minutes before measurement.
- Avoid taking a measurement for 30 minutes after physical exercise, bathing or eating.
- For greater measurement accuracy, at least 2 minutes should pass between two measurements.
- Avoid direct sunlight or other heat sources, or direct air conditioner flow on the forehead while taking temperature, since these factors might interfere with the measurement, influencing results.


When pointed at a hot object, the thermometer does not measure the actual temperature of the object, since IR radiation depends on the characteristics of the body being measured and the device in question is obviously calibrated on the coefficient of radiation of human skin (which is very different from that of any other material). To avoid erroneous readings, minimize the time from thermometer "READY" to scanning the forehead.

- Do not expose the thermometer to extreme temperatures (lower or higher than the prescribed storage temperatures -20°C + 55°C), dust or shock; otherwise the sensor could be irreparably damaged.
- Do not splash the thermometer with water.

Features of the Thermo-Precision thermometer

1. Optical module
2. LCD display
3. Measurement button
4. Memory button
5. On/off ring
6. Battery compartment

FIGURE 1

Measurement

1. To switch on the thermometer, turn the ring clockwise, until completely open, i.e. the ring notch should coincide with the notch of the thermometer (FIGURES 2 and 3). It is important that the ring is turned completely, in order to obtain the full opening of the thermometer "objective", causing the blue leds to switch on and their light to converge in a single blue ray that indicates the precise point where the temperature is read. The blue ray is very useful to ensure the correct positioning of the probe during the reading. To start the reading, wait for the sign "READY" to appear. (FIGURE 4).


FIGURE 2


FIGURE 3

FIGURE 4

2. Hold the thermometer at a distance not exceeding 10 cm from the forehead. Press the "Measurement" button (FIGURE 1) and scan the whole temple area. The best position is at the left and above the left eyebrow or at the right and above the right eyebrow, just underneath the hairline.

While scanning you will hear short beeps which indicate that measurement is taking place.

3. The measurement is complete only after 4 beeps have sounded. At this point you can read the temperature on the display. The temperature appears now on the display.
4. The thermometer switches automatically off 1 minute after measurement has been completed. To switch off the guiding blue ray, rotate the ring (if the thermometer "objective" is not closed completely, the ray remains active, causing the batteries that power it to become rapidly exhausted).

IMPORTANT: Any distance from 0.5 cm up to 10 cm, with the thermometer correctly positioned, provides correct measurements.

IMPORTANT: For repeated measurements always use the same measuring area and side of the forehead.

Memory function

The memory function saves 8 measurements for 64 hour after each measurement.

The last reading before the thermometer switches automatically off, is saved in the memory. To use the memory function, turn the ring (FIGURE 2) and press and release the Memory button on the front side of the thermometer (FIGURE 1) while the thermometer is on. The LCD display alternately displays the following information four times:

- Previous measurement.
- How many hours before the measurement was taken (NOTE: only whole hours are indicated, until the first full hour has elapsed, the value appearing on the display is 0).

Pressing and releasing again the Memory button, the measurement prior to the last is displayed.

You can read all your saved readings by continuously pressing the memory button.

To cancel all saved measurements, press the Memory button while turning the ring anticlockwise; in this way the thermometer is switched off and all data will be cancelled.

DEGREES CELSIUS (°C) / DEGREES FAHRENHEIT (°F) SWITCH

To change the measurement unit from degrees Celsius to Fahrenheit, press the Memory button while the thermometer is switched off.

The display shows °C or °F, according to the selected mode, for 5 seconds. To switch from °C to °F and vice versa, press again the Memory button while the display is on.

BATTERY REPLACEMENT

The forehead infrared Thermo-Precision thermometer is supplied with a 3V lithium battery, already installed. Battery shelf life is five years. Battery's working life is approximately 12000 measurements.

Besides, the thermometer contains two 1.5V AAA batteries powering the blue ray.

If the "Low Battery" symbol is displayed

or if the exhausted battery symbol appears

On the thermometer's LCD display or the thermometer does not function at all, immediately replace the battery as follows:

- Pull the battery cover straight back and upwards to remove it.
- Remove the old battery from the compartment.
- Insert the new battery with the positive (+) side up and the negative (-) side pointed down until it snaps into place with an audible "click".
- Replace battery cover.

Note: The "Low Battery" symbol is intended for 3V lithium (operating) battery only. The AAA batteries should be replaced in case of blue light just not appears.



CAUTION: replace only with CR 2032 batteries (for example: Energizer, Varta, Panasonic).

CAUTION: The use of other types of batteries may present a risk of fire or explosion.

CAUTION: Do not recharge, disassemble, heat above 100 °C (212°F), or incinerate.

If the thermometer does not project the blue ray, replace the two AAA batteries as follows:

- remove the battery cover;
- remove the old batteries;
- insert the new batteries with the positive (+) and the negative (-) poles arranged as indicated by the icons in the compartment.

Caution: replace only with AAA batteries (e.g.: Duracell, Energizer).

The use of other types of batteries may cause explosions or fire.

Caution: Battery may explode if mishandled. Do not recharge, disassemble, heat above 100 °C (212°F), or incinerate.

CARE AND MAINTENANCE

The thermometer has an infrared sensor inside the optical module. To avoid damage to the sensor, when the thermometer is not used, close the thermometer by turning the ring anticlockwise. If dust, water or other traces of dirt accumulate inside the optical module, or if it shows finger smudges on its surface, gently clean it with a soft cloth and allow to dry before closing it again. It is advised to occasionally clean the external part of the optical module (without touching the sensor) with a 70% alcohol solution.

TROUBLESHOOTING

Display	Meaning	Action
	There is a general thermometer function problem. Note: the small digits from 0 to 9 in the circle are intended for internal reference of manufacturer, but not for a customer	Reset the thermometer by pressing the "ON" and "M" buttons simultaneously and keeping them pressed for 3 seconds. If the thermometer still does not function, contact technical services.
	The operating temperature is lower than 16°C (60.8°F).	Place the thermometer in a room where the temperature is between 16°C and 36°C (60.8°F - 96.8°F). Wait 30 minutes before taking another temperature reading.
	The operating temperature is higher than 36°C (96.8°F).	Place the thermometer in a room where the temperature is between 16°C and 36°C (60.8°F - 96.8°F). Wait 30 minutes before taking another temperature reading.
	The temperature measured is instable	Wait 2 minutes and repeat the measurement strictly as described in this manual above.
	The temperature measured is outside the measurement range 34.4°C - 42.6°C (93.4°F - 108.7°F).	Wait 2 minutes, then repeat the measurement and wait for the 4 beeps before moving the thermometer.
	The battery is low. Only a few more measurements can be taken.	Replace the old battery with a new one.
	The battery is dead. No more measurements can be taken.	Replace the old battery with a new one immediately.
The thermometer does not function or does not respond. The display does not go on or remains on steady.	The thermometer needs to be reset or the battery replaced.	Reset the thermometer by pressing the "ON" and "M" buttons simultaneously and keeping them pressed for 3 seconds or replace the old battery with a new one. If the thermometer still does not function, contact technical services.

Note: In case of any error Buzzer sounds one beep duration 1s

TECHNICAL SPECIFICATIONS

Temperature measurement range: 34.5°C-42.6 °C (94.1 °F-108.7°F)

Operating temperature range: 16.0°C-36.0°C (60.8 °F-96.8°F)

Laboratory accuracy: $\pm 0.2^{\circ}\text{C}$ ($\pm 36^{\circ}\text{F}$)

Batteries:

one 3V Lithium Battery CR 2032 (for measurement functions)

two 1.5V type AAA batteries (only for the blue ray)

Clinical accuracy:

Clinical Bias: $< 0.03^{\circ}\text{C}$

Clinical Repeatability: $\leq \pm 0,3^{\circ}\text{C}$

Accuracy according to ASTM E-1965-98 & EN 12470-5 standards.

Clinical accuracy characteristics and procedures are available on request.

Classification according to safety regulation EN60601-1:

1. Internally powered equipment.
2. Type BF.
3. Continuous Operation

Caution: keep this instruction in a safe place and read them before use.

The thermometer should be stored at room temperature in a dry place, away from heat sources and direct sunlight. Strong electromagnetic fields might affect its correct operation. Keep out of reach of children. Never leave a child unattended while taking its temperature.

In case of a rise in temperature, consult a physician. Conforming to directive 93/42 EEC.

2 years warranty

Artsana's Thermo-Precision thermometer is guaranteed to be free of defects due to the use of faulty materials or manufacturing faults for a period of two years from the date of purchase. This warranty applies only if the device is used at normal domestic conditions in conformity to the instruction manual.

In no case Artsana is liable for particular, direct, indirect or consequent damage deriving from or related to the purchase or use of this product, or for higher costs than the original price of the product.

The warranty does not cover any damage due to the battery or accidental damage caused by or resulting from improper use of the device. Further warranties may vary from state to state.

ELECTROMAGNETIC COMPATIBILITY

This device is manufactured in conformity to the regulations in force and complies with all protection requirements in terms of electromagnetic compatibility; nevertheless, special precautions should be taken concerning EMC and the device should be installed and operated in conformity to the EMC information contained in this document.

Portable and mobile radiocommunication devices may influence the operation of this electromedical device. In case of electromagnetic disturbance, it is recommended to repeat the measurement of body temperature, if the readings appear dubious.

A guide to the manufacturer's statement – electromagnetic emissions

The digital electronic thermometer Thermo-Precision has been manufactured to operate in the following electromagnetic environment. The customer or the user of the thermometer should make sure that it is used in this environment.

Emission test	Conformity	Electromagnetic environment - guide
RF emissions	Group 1	The Thermo-Precision thermometer uses RF energy only for its internal functioning. Therefore its RF emissions are very low and presumably do not cause any interference with neighbouring electronic devices.
RF emissions	Class B	The Thermo-Precision thermometer is suitable for use in all rooms, including domestic rooms and rooms directly connected to the public low voltage mains supplying civil buildings.

0123

REF 00071691000000

Manufactured by Sola Electro –Optics Co. Ltd. Shanghai China

EC representative: Obelis s.a. Av. de Tervuren, 34 Bte 44 B-1040 Brussels, Belgium

Distributed by: Artsana S.p.A. - Grandate (CO) Italy www.chicco.com

This product complies with EU Directive 2002/96/EC.

The crossed bin symbol on the appliance indicates that the product, at the end of its life, must be disposed of separately from domestic waste, either by taking it to a separate waste disposal site for electric and electronic appliances or by returning it to your dealer when you buy another similar appliance. The user is responsible for taking the appliance to a special waste disposal site at the end of its life. If the disused appliance is collected correctly as separate waste, it can be recycled, treated and disposed of ecologically; this avoids a negative impact on both the environment and health, and contributes towards the recycling of the product's materials. For further information regarding the waste disposal services available, contact your local waste disposal agency or the shop where you bought the appliance.

THIS PRODUCT COMPLIES WITH EC DIRECTIVE 2006/66/EC

The crossed bin symbol on the batteries indicates that, at the end of their life, they must be disposed of separately from domestic waste, either by taking them to a separate waste disposal site for batteries or by returning them to your dealer when you buy similar rechargeable or non-rechargeable batteries. The chemical symbols Hg, Cd, Pb, printed under the crossed bin symbol, indicate the type of substance contained in the batteries: Hg=Mercury, Cd=Cadmium, Pb=Lead. The user is responsible for taking the batteries to a special waste disposal site at the end of their life, so that they can be treated and recycled. If the spent batteries are collected correctly as separate waste, they can be recycled, treated and disposed of ecologically; this avoids a negative impact on both the environment and human health, and contributes towards the recycling of the batteries' substances. Non-compliance with the norms on battery disposal damages the environment and human health. For further information regarding the waste disposal services available, contact your local waste disposal agency or the shop where you bought the batteries.

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Artsana S.p.A. - Via Saldarini Catelli, 1
SOLA ELECTRO-OPTICS.co.Ltd